

## SOP\_001\_NU\_1\_2\_Top\_Down\_Standard\_v2D\_CJD\_LF

Caroline J. DeHart and Luca Fornelli

## ❖ Reagent and Materials List

Item	Part Number	Vendor
Ubiquitin	U6253	Sigma Aldrich
Trypsinogen	T1143	Sigma Aldrich
Myoglobin	M5696	Sigma Aldrich
Carbonic Anhydrase	C2624	Sigma Aldrich
Optima Grade Water	W6-4	Fisher Scientific
Optima Grade Acetonitrile	A955-4	Fisher Scientific
MS-Grade Formic Acid	PI-28905	Fisher Scientific
1.5 mL Protein LoBind Microcentrifuge Tubes	13-698-794	Fisher Scientific

## ❖ Important Notes

- ◆ Use 1.5 mL Eppendorf LoBind microcentrifuge tubes for protein stock preparation, top-down (TD) standard preparation, and long term aliquot storage. In our experience, these tubes have shown the lowest degree of plasticizer leaching and/or protein binding during use and storage.
- ◆ Approximate final protein amounts (loaded on-column): 0.1 pmol ubiquitin, 0.5 pmol trypsinogen, 1 pmol myoglobin, and 0.6 pmol carbonic anhydrase. Superoxide dismutase (SOD1) is present as a contaminant in carbonic anhydrase.
- ◆ A TD standard prepared in this way should be stable for up to three days at 4 °C (before significant protein oxidation becomes evident).

## ❖ Recipe

- ◆ Prepare 2 mg/mL stocks of each protein standard in Optima H<sub>2</sub>O. (Aliquots can be stored at -80 °C.)
- ◆ Prepare the following (volumes shown from respective stock solutions):

Protein	Volume (μL)	Stock Concentration (pmol/μL)	Amount Loaded on Column (pmol, 1X)
Carbonic Anhydrase	40	25.7	0.64
Myoglobin	40	43.9	1.09
Trypsinogen	25	19.6	0.49
Ubiquitin	2.5	5.5	0.14
<b>Total</b>	<b>107.5</b>		

- ◆ Divide final mixture into 2.5 uL aliquots and store at -80 °C.

## ❖ Preparation

- ◆ Dilute one aliquot of TD STD in 100x vol. of Buffer A (95% Optima H<sub>2</sub>O, 5 % Optima Acetonitrile, 0.2% MS-grade formic acid), where 1x vol. is the intended injection volume (e.g. **600 µL Buffer A** for an intended injection volume of **6 µL**). This will ensure that the correct amount of each TD standard protein is present in each injection.
- ◆ Mix thoroughly by pipetting, then transfer to a clean autosampler vial. The standard is now ready for use.


## ❖ Column Parameters

- ◆ **Option 1:** Self-packed PLRP-S columns
  - Packing Material: PLRP-S resin, 1000 Å pore size, 5 µm particle size (obtained from Agilent Technologies)
  - Trap column: 2 cm bed length, 150 µm I.D.
  - Analytical column: 15 cm bed length, 75 µm I.D.
  - Nanospray Emitter: 15 µm PicoTip emitter, packed with 2 mm PLRP-S resin (P/N FS360-50-15-N-20-C12, New Objective)
- ◆ **Option 2:** Thermo Dionex Monolithic RP-4H columns
  - Trap column: 0.5 cm bed length, 150 µm I.D. Thermo Dionex PepSwift (P/N 164558)
  - Analytical column: 50 cm bed length, 100 µm I.D. Thermo Dionex ProSwift Monolithic RP-4H (P/N 164921)

## ❖ LC Parameters

- ◆ **Solvent A:** 95% Optima H<sub>2</sub>O, 5% Optima Acetonitrile, 0.2% MS-grade formic acid
- ◆ **Solvent B:** 5% Optima H<sub>2</sub>O, 95% Optima Acetonitrile, 0.2% MS-grade formic acid
- ◆ **Option 1:** Self-packed PLRP-S columns
  - Trapping configuration: 3 µL/min flow rate (10 min. trap cycle, 45 °C)
  - Analytical configuration: 0.3 µL/min flow rate (48 min. analytical gradient, 45 °C)
  - Gradient Parameters:


Time (min.)	% B	Curve
0.0	5.0	
10.0	5.0	0% in 10 min.
13.0	15.0	10% in 3 min.
37.0	50.0	35% in 24 min.
40.0	95.0	45% in 3 min.
42.0	95.0	0% in 2 min.
45.0	5.0	90% in 3 min.
58.0	5.0	0% in 13 min.



◆ **Option 2:** Thermo Dionex Monolithic RP-4H columns

- Trapping configuration: 10 µL/min flow rate (3 min. trap cycle, 35 °C)
- Analytical configuration: 1 µL/min flow rate (55 min. analytical gradient, 35 °C)
- Gradient Parameters:

Time (min.)	% B	Curve
0.0	5.0	
3.0	5.0	0% in 3 min.
33.0	55.0	50% in 30 min.
35.0	95.0	40% in 2 min.
38.0	95.0	0% in 3 min.
39.0	5.0	90% in 1 min.
53.0	5.0	0% in 14 min.



❖ **MS Parameters**

◆ **Instrument Tuning and Method Parameters:**

(All in positive and profile mode, with 15.0 V source CID and “Protein Mode” on.)


<b>Scan Event 1: FTMS1</b>	Scan Range ( <i>m/z</i> )	500.00 – 2000.00
(120k RP)	Microscans	4
Full Scan	Max Inject Time (ms)	50.00
Normal mass range	MS1 AGC Target	1.00e +06
<b>Scan Event 2: FTMS2</b>	Activation type	HCD
Scan range	Define <i>m/z</i> range	400-2000
(60k RP)	Default charge state	10
Quadrupole isolation ON	Isolation Width ( <i>m/z</i> )	3.0
Top 2, dd	Normalized Collision Energy	22.0
	Microscans	4
	Max Inject Time (ms)	800.00
	MS2 AGC Target	1.00e+06

◆ **Dynamic Exclusion Settings (MS2):**

Underfill Ratio	2 %
Exclude Isolation (Charge states 1-5)	On
Exclusion Duration (s)	100
Intensity Threshold	2.50e+04


❖ Data Interpretation and Analysis

◆ Example Chromatograms:




**Example Chromatograms:** Typical base peak chromatogram of the NRTDP TD standard on the PLRP-S (left) and RP-4H (right) columns described above, showing four separate eluted protein peaks (1. Ubiquitin, 2. Trypsinogen, 3. Myoglobin, 4. Carbonic Anhydrase). Superoxide dismutase, a characteristic contaminant of carbonic anhydrase, is indicated in the RP-4H chromatogram by an asterisk. The elution order and relative height ratio of all other protein peaks should remain consistent. The examples shown were obtained on a QE-HF, using the LC and MS parameters described above.

◆ Example FT MS1 spectra:



**Example FT MS1 spectra:** Averaged FTMS1 spectra for each of the four peaks in the above PLRP-S and RP-4H chromatograms, showing the characteristic isotopic peak distributions for each protein.


◆ Example FT MS2 spectra:



**Ubiquitin**  
(UniProt ID P0CG53)

```
N M[Q]I[F]V[K]T L[T]G K[T]I[T]L[L]E[V]E[P]S[D]T[I]E[N] 25
26 [V]K A[K]I[Q]D[K]E[G]I[P]L[D]Q Q R L I F A G K[Q]L 50
51 [E]D[G]R[T]L S D[Y]N[I]Q[K]E[S]T[L]H L V L R L R G 75
76 G C
```


P-score: 4 E-102, PCS: 1226, 75 matching fragments



**Myoglobin**  
(UniProt ID P68082)

```
N G L S D G E[W]Q[Q]V[L]N V W G K V E A D I A G H G 25
26 Q E V L I R L F T G H P E T L E K F D K F K H L K 50
51 T E A E M K A S E D[L]K K H[G]T V V L T A L G G I 75
76 L K K K G H H E A E L K P L A Q S H[A]T K[H]K I[P] 100
101 I[K]Y[L]E F[I]S D A I[I]H[V]L H S K[H]P G D[F]G A 125
126 [D]A[Q]G A[M]T K A[L]E[L]F R N[D]I[A]A[K]Y[K]E[L]G 150
151 F Q G C
```

P-score: 2 E-33, PCS: 335, 42 matching fragments



**Carbonic Anhydrase**  
(UniProt ID P00921)

```
N S H H[W]G[Y]G[K]H N[G]P E[H]W H K[D]F[P]I[A]N G E 25
26 R[Q]S[P]V[D]I[D]T K A V[V]Q D[P]A[L]K P L A L V Y 50
51 G E A T S R R M V N N G H S F N[V]E Y D D S Q D K 75
76 A V L K D G P L T G T Y R L V Q F H F H W G S S D 100
101 D Q G S E H T V D R K K Y A A E L H L V H W N T K 125
126 Y G D F G T A A Q Q P D G L A V V G V F L K V G D 150
151 A N P A L Q K V L D A L D S I K T K G K S T D F P 175
176 N F D P G S L L P N V L D Y W T Y[P]G S L T T[P]P 200
201 L L E S V[T]W I V[L]K[E]P I[S]V S S[Q]Q[M]L[K]F R 225
226 T L N F N A[E]G[E]P[E]L[L]M[L]A[N]W[R]P A Q P L K 250
251 N R Q V R[G]F P K C
```

P-score: 6 E-67, PCS: 766, 51 matching fragments

**Example FT MS2 spectra:** Single-scan fragmentation spectra for ubiquitin (**top**), myoglobin (**middle**), and carbonic anhydrase (**bottom**) from the RP-4H dataset shown above. The fragment ion masses from each of the above spectra were deconvoluted using the Xtract algorithm (Thermo) and searched against the respective protein sequences using ProSight Lite. Typical P-scores for ubiquitin should be below E -50, while P-scores for myoglobin and carbonic anhydrase should be below E -25.

## ◆ Data Analysis Methods:

- **ProSight Lite:** The software is available for free download at <http://prosightlite.northwestern.edu/>. A detailed protocol for the analysis of the NRTDP Top-Down Standard with Xtract and ProSight Lite can be found at [https://link.springer.com/content/pdf/10.1007%2F978-1-4939-6783-4\\_18.pdf](https://link.springer.com/content/pdf/10.1007%2F978-1-4939-6783-4_18.pdf)
- **ProSight PC 4.0:** A “Standards” search database for high-throughput data analysis of the NRTDP Top-Down Standard with ProSight PC 4.0 is available for download here: <http://proteinaceous.net/database-warehouse/>
- **NRTDP TDPortal:** A custom workflow for high-throughput analysis of the NRTDP Top Down Standard is available on the TDPortal Quest-based, high-performance computing environment available through NRTDP and Northwestern University. User accounts can be requested at <http://nrtdp.northwestern.edu/tdportal-request/>. A detailed protocol for data analysis on TDPortal by external users (**NRTDP SOP\_004**) can be found at <http://nrtdp.northwestern.edu/wp-content/uploads/2017/01/ExternalUserJan10.pdf>

## ❖ Longitudinal Data Tracking

The NRTDP recommends including the following metrics into longitudinal tracking of LC and MS performance:

- ◆ **Peak Intensity:** IT and FT MS1 peak intensity of ubiquitin, myoglobin, and carbonic anhydrase
- ◆ **Peak Area (log):** Chromatographic peak area (and retention time) of ubiquitin, myoglobin, and carbonic anhydrase
- ◆ **FWHM:** Full width at half maximum of ubiquitin, myoglobin, and carbonic anhydrase peaks
- ◆ **Injection Time:** MS1 and MS2 injection times for ubiquitin, myoglobin, and carbonic anhydrase
- ◆ **P-score (-log):** ubiquitin, myoglobin, and carbonic anhydrase obtained by low- or high-throughput data analysis

Paying close attention to these parameters over time can help identify LC or MS issues before they become significant, thus reducing loss of important sample data.

The NRTDP further recommends running at least three injections of the Top Down Standard before and after running experimental samples, as well as at least one injection of Top Down Standard every twenty-four hours. Evaluation of these standards on the fly will not only help detect LC or MS issues, but also provide confirmation that optimal performance is maintained when consistency and reproducibility are crucial.

**Acknowledgements:** This protocol was developed with support from the National Institute of General Medical Sciences **P41GM108569** for the National Resource for Translational and Developmental Proteomics (NRTDP) based at Northwestern University.